

B/S BRACE SYSTEM

The "classic" for problem-free use

The B/S Brace is a thin, easy-to-use plastic adhesive bar for cosmetic treatment of sensitive toenails. It is glued over the entire nail surface, optimally distributing the necessary spring forces to provide relief for the affected tissue. Both the optimal fit and firm adhesion allow the patient a high level of wearing comfort.



1. Thorough preparation



2. Careful degreasing



3. Exact fitting from the start



4. Exact fitting to the end



5. Grinding the ends of the brace



6. Application of the B/S activator



7. Application of the B/S Basic adhesive



8. Correct centre fixing



9. Pressing the B/S Brace into place on the right



10. Pressing the B/S Brace into place on the left



11. Aligning the brace transitions



12. Individual tension correction



13. Careful final cleaning



14. Final sealing



15. Satisfying results